


Falke with Flow Stand. O2 DIN 9/16



Manufacturer: Greggersen

Product 910350

Number:

Description

Falke with Flow, Standard

Compact pressure regulator with contents gauge;
incrementally adjustable flow rate.

Gas Type: Oxygen / O2

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 9

Outlet: 9/16-18 UNF (tubing nipple) bottom

Flow Rate: 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type O2-Oxygen