


Forano, Base Part, VS100M A rh, DIN, BS, VAC



Manufacturer: Greggersen

Product 903212

Number:

Description

Terminal Unit Forano Base Part DIN/AGA/NF/BS

Gas-specific base part for installation in

Trilux VS100 media supply duct.

Base part with spacer washer and mounting bracket;

brazed copper pipe, integrated non-return and maintenance valve.

Scope of delivery:

- gas-specific base part with Cu pipe
- spacer washer (2 pcs)
- mounting bracket

Technical data:

Gas type: Vacuum / VAC

Standard AGA: UTL (AGSS)

Installation type: Trilux VS100 M (A)

Inlet: Cu pipe 8x1 mm bent, right

Properties

Gas Type

VAC-Vacuum
