

Kolibri, DIN Probe O2 Standard ISO32



Manufacturer: Greggersen

Product 902160

Number:

Description

Kolibri, O₂, 0...15 l/min

- Standard/PLUG-IN UNIT -

Medical compact flowmeter according to DIN EN ISO 15002,
for connection to a terminal unit, with handwheel
step-adjustable flow rate, outlet 9/16-18 UNF
and tubing nipple

Handwheel in ISO32 - white for O₂

Gas type: Oxygen / O₂

Connection: Probe according to DIN 13260 Part 2

Total length: 88 mm (probe front part to centre of tubing nipple outlet)

Outlet: 9/16-18 UNF / tubing nipple

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
