


Angled Probe O2 Hose Barb 7mm, DIN, ISO V4A



Manufacturer: Greggersen

Product 505150

Number:

Description

Angled Probe for Oxygen

with 6 mm Hose Barb and Non-Return Valve

Probe Front Part made of 1.4404 (Stainless Steel)

Gas Type: Oxygen / O2

Connection: Probe according to DIN 13260 Part 2

Outlet: 7 mm Hose Barb

Colour Coding: according to ISO 9170-1 (white)

Properties

Gas Type

O2-Oxygen