

Double Kolibri, NIST, O2 Standard



Manufacturer: Greggersen

Product 902678

Number:

Description

O2 - Double Kolibri, 0...15 l/min Standard

RAIL-MOUNTED UNIT NIST

Medical Compact Flowmeter according to DIN EN ISO 15002;

flow rate is adjusted via a handwheel in 10 flow settings; gas inlet NIST integrated into rail clamp; outlet 9/16-18 UNF and tubing nipple

Gas type: Oxygen O2

Connection: NIST according to DIN EN ISO 18082

Outlet: 2 x 9/16-18 UNF / tubing nipple

Flow Rate: 2 x 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	O2-Oxygen
-----------------	-----------