

Kolibri, BS Probe, O2 Standard



Manufacturer: Greggersen

Product 901440

Number:

Description

Kolibri, O2, Standard, BS

PLUG-IN UNIT

Medical compact flowmeter according to
DIN EN ISO 15002, for connection
to a terminal unit, with handwheel
incrementally adjustable flow rate,
outlet 9/16-18 UNF and tubing nipple

Gas type: Oxygen

Connection: Probe type BRITISH STANDARD

Outlet: 9/16-18 UNF / tubing nipple

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
