

Habicht Rail, O2, NIST, Kolibri Antimicrobial



Manufacturer: Greggersen

Product 500092

Number:

Description

Emergency Supply Unit

Suction and Insufflation Unit, for Oxygen

RAIL-MOUNTED UNIT NIST

ANTIMAGNETIC

Combi-unit for attachment to an equipment rail;
gas supply via a pressure regulator or
a gas terminal unit; ejector regulator with
vacuum gauge, fine adjustment valve and silencer
with filter holder; Kolibri compact flowmeter with
10 flow settings

Gas type: Oxygen / O2

Connection: NIST in accordance with DIN EN ISO 18082

Flowmeter outlet: 9/16-18 UNF with Tubing Nipple

Ejector outlet: Vacuum Tubing Nipple

consisting of:

- Rail clamp with NIST pressure inlet
- Distribution block, spacer
- Kolibri Standard 0...15 l/min
- Ejector regulator, 0 to -0.9 bar