

Kolibri, DIN Probe, O2 Standard antim.



Manufacturer: Greggersen

Product 500231

Number:

Description

Kolibri, O₂, 0...15 l/min, antimagnetic

- Standard/Plug-in Unit -

Medical Compact Flowmeter according to DIN EN ISO 15002,

for coupling into a terminal unit, with handwheel

stepwise adjustable flow rate,

outlet 9/16-18 UNF and tubing nipple

Version: Antimagnetic

Gas Type: Oxygen

Connection: Probe according to DIN 13260 Part 2

Overall Length: 88 mm (probe front part to centre of tubing nipple outlet)

Outlet: 9/16-18 UNF / Tubing Nipple

Flow Rate: 0...15 l/min

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
