

Kolibri 2 Probe MEDAP O2 Standard



Manufacturer: Greggersen

Product 906450

Number:

Description

KOLIBRI 2, O₂, 0...15 l/min

- Standard MEDAP PLUG-IN UNIT -

Flow settings:

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selector device

Medical compact flowmeter for connection

to a terminal unit, with handwheel - incrementally

adjustable flow rate, outlet 9/16-18 UNF thread

Adjustment device ergonomically optimised and readability

of the set flow setting possible from various viewing directions

and angles.

Frontal operating concept -> increased safety when
setting the correct flow rate

Gas type: Oxygen / O₂

Connection: Probe suitable for MEDAP

Outlet: 9/16-18 UNF

Flow rate: 0...15 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:

900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	O ₂ -Oxygen
Inlet Standard	MEDAP
Variant	Standard