

AGSS Forano AIR UP/HW Type1L



Manufacturer: Greggersen

Product 902073

Number:

Description

AGSS FORANO Type 1

for flush and cavity wall installation

according to DIN EN ISO 9170-2, with ejector,
exhaust air return, with visual indicator
for operational check, in stainless steel housing
(146 x 135 x 64 mm), and front plate made of
stainless steel (174 x 174 mm)

Flow rate: min. 50 l/min at 500 kPa (5 bar)

Gas type: Compressed air / AIR

Inlet: Cu pipe Ø 8 mm

Outlet: Exhaust air Cu pipe Ø 15 mm

Properties

Gas Type

AIR Compressed Air
