

Kolibri 2 Double Rail NIST O2 Standard



Manufacturer: Greggersen

Product 906428

Number:

Description

KOLIBRI 2, O₂, 0...15 l/min

- Standard - Double - RAIL-MOUNTED UNIT -

Flow Settings (2x):

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Double Compact Flowmeter KOLIBRI 2

with fixed orifices and selection device

Medical compact flowmeter for mounting on an
equipment rail 25x10mm; gas inlet NIST integrated in the
rail clamp; with handwheel - incrementally
adjustable flow rate, outlet 9/16-18 UNF thread

Adjustment device ergonomically optimised and readability
of the set flow setting from various viewing directions

and angles possible.

Frontal operating concept -> greater safety when
setting the correct flow rate

Gas type: Oxygen / O2
Connection: NIST in accordance with DIN EN ISO 18082
Outlet: 9/16-18 UNF
Flow rate: 0...15 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:
2x 900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	O2-Oxygen
Inlet Standard	NIST
Variant	Standard