

Kolibri, DIN Probe, CO2 Standard



Manufacturer: Greggersen

Product 901419

Number:

Description

Kolibri, CO2, 0...15 l/min

Standard/PLUG-IN UNIT

Medical compact flowmeter according to DIN EN ISO 15002,
for coupling into a terminal unit, with handwheel
step-adjustable flow rate,
outlet 9/16-18 UNF and tubing nipple

Gas type: Carbon dioxide / CO2

Connection: Probe according to DIN 13260 Part 2

Overall length: 88 mm (probe front part to centre of tubing nipple outlet)

Outlet: 9/16-18 UNF / tubing nipple

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	CO2-Carbon Dioxide
----------	--------------------
