

Double Flowmeter, DIN, O₂, 0-15l/min



Manufacturer: Greggersen

Product 900745

Number:

Description

Double Flowmeter

-PLUG-IN UNIT-

Medical double flowmeter,
pressure-compensated according to DIN EN ISO 15002,
with fine control valve for dosing the flow,
for coupling into a terminal unit

Gas type: Oxygen / O₂

Connection: Probe according to DIN 13260 Part 2

Outlet: 2 x 9/16-18 UNF

Flow rate: 0...15 l/min each, continuously adjustable

Inlet pressure: 450 kPa +/- 50 kPa

Properties

Gas Type

O2-Oxygen
