

Kolibri 2 Probe DIN O2 Standard+ Long



Manufacturer: Greggersen

Product 906406

Number:

Description

KOLIBRI 2, O₂, 0...30 l/min

- Standard+ DIN PLUG-IN UNIT -

Long version

Flow settings:

0 / 2 / 3 / 4 / 6 / 8 / 12 / 15 / 20 / 25 / 30 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices with selector mechanism

Total length "long"

Probe front part to centre of flow outlet

(Distance from terminal unit to tubing nipple):

75 mm (intermediate piece 60mm)

Medical compact flowmeter for plugging

into a terminal unit, with handwheel - incrementally

adjustable flow rate, outlet 9/16-18 UNF thread

Adjustment mechanism ergonomically optimised and readability of the selected flow setting possible from various viewing directions and angles.

Frontal operating concept -> greater safety when setting the correct flow rate

Gas type: Oxygen / O2
Connection: Probe according to DIN 13260 Part 2
Outlet: 9/16-18 UNF
Flow rate: 0...30 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:
900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	O2-Oxygen
Inlet Standard	DIN
Variant	Standard+ Long