


Falke without Flow easy O2 DIN 3/8"



Manufacturer: Greggersen

Product 910372

Number:

Description

Falke "easy", 450 kPa, G3/8",

Compact pressure regulator with high-pressure cylinder connection to DIN 477, contents gauge, free flow at a pressure of 450 kPa.

Gas type: Oxygen / O₂

Connection: High-pressure cylinder connection to DIN 477, No. 9

Outlet: G3/8" with non-return valve, bottom

Output: 450 kPa

Properties

Gas Type	O ₂ -Oxygen
-----------------	------------------------