

Flowmeter, BS, O₂, 0-15l/min



Manufacturer: Greggersen

Product 901501

Number:

Description

Flowmeter, O₂, 0...15 l/min, BS

- PLUG-IN UNIT -

Medical flowmeter, pressure-compensated according to DIN EN ISO 15002, with fine adjustment valve for dosing the flow, for connection to a terminal unit

Gas type: Oxygen

Connection: Probe type BRITISH STANDARD

Outlet: 9/16-18 UNF

Flow rate: 0...15 l/min continuously adjustable

Inlet pressure: 450 kPa +/- 50 kPa

Properties

Gas Type

O2-Oxygen
