


Ventus 22mm Valve Insert AIR Sensor16 b



Manufacturer: Greggersen

Product 900779

Number:

Description

Valve Insert for Ventus 22mm

Gas type: AIR

Design complies with DIN EN ISO 7396-1

Valve assembly consisting of:

- Shut-off valve DN 20, G3/4"
- Copper pipe 22mm inlet / outlet
- Pressure gauge 0...16 bar
- Pressure sensor 0... 16 bar
- NIST emergency inlet point
- Physical separation point
- Gas type labelling

Properties

Fulfilled Standard	DIN EN ISO 7396-1
Gas Type	AIR Compressed Air
NIST Emergency Inlet Point	1 NIST Emergency Inlet Point