


Kolibri, MEDAP Probe, AIR Standard



Manufacturer: Greggersen

Product 901584

Number:

Description

Kolibri, Standard

Medical compact flowmeter according to DIN EN ISO 15002,
for coupling into a terminal unit, with handwheel
stepwise adjustable flow rate,
outlet 9/16-18 UNF and tubing nipple

PLUG-IN UNIT

Gas type: Compressed Air / AIR

Connection: MEDAP Probe

Outlet: 9/16-18 UNF / Tubing Nipple

Flow rate: 0...15 l/min

Properties

Gas Type AIR Compressed Air