

Forano, Base Part, VS100M A re,DIN,O2



Manufacturer: Greggersen

Product 903210

Number:

Description

Terminal Unit Forano Base Part (Base Part) DIN/AGA/NF

Gas-specific base part for installation in

Trilux VS100 media supply unit.

Base part with spacer washer and mounting bracket;

brazed copper pipe, integrated non-return and maintenance valve.

Scope of delivery:

- Gas-specific base part with Cu pipe
- Spacer washer (2 pcs.)
- Mounting bracket

Technical data:

Gas type: Oxygen / O₂

Installation type: Trilux VS100 M (A)

Inlet: Cu pipe 8x1 mm bent, right

Properties

Gas Type

O2-Oxygen
