


Forano, Front Part, DIN, AIR-800, ISO



Manufacturer: Greggersen

Product 905109

Number:

Description

Terminal Unit FORANO - Front Part DIN

ISO 32 -

Gas type: Compressed Air-800 / AIR-800

Inscription: medical air 800 / AIR-800

Gas-specific socket assembly for mounting
on FORANO base part.

Socket assembly according to DIN 13260 Part 2
with parking and operating position. Release via release ring.

Colour coding:

according to DIN EN ISO 9170-1 (ISO32 - black white)

Release ring inscription: printed

Language: English

Cartridge incl. sealing elements

Technical Data:

Gas type: Compressed Air 800 / AIR-800

Inscription: medical air 800 / AIR-800

Colour coding: according to ISO 9170-1 (black-white)

Front diameter: Ø 42mm

Scope of Delivery:

- Gas-specific socket assembly
- Release ring incl. printed inscription
- Cartridge with sealing elements

Properties

Gas Type

AIR-800-Compressed Air 8bar
