

Hose, AIR, 3m, DIN - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 900362

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow restrictor,
hose and NIST device connection

Specification:

Gas type: Compressed air / AIR

Inlet: Probe according to DIN 13260 Part 2

Hose length: 3.0 m

Colour coding: Colour-Neutral

Outlet: gas-specific NIST screw connection

Properties

Gas Type AIR Compressed Air