


Medical Hose 6x12 Black/White ISO AIR



Manufacturer: Greggersen

Product 900342

Number:

Description

Connecting hose for med. gases

Braided hose

Inner diameter: 6.7 +/- 0.15 mm

Wall thickness: 3.0 +/- 0.2 mm

Gas type: Compressed Air

Colour coding: ISO AIR

for med. gases, black/white striped

Text for outer hose (white):

AIR...125mm...AIR...125mm...AIR

PU: Roll of 100 metres

Properties

Gas Type

AIR Compressed Air