


Falke with Flow Standard O2 DIN antimag +1DIN



Manufacturer: Greggersen

Product 500158

Number:

Description

Falke with Flow, Standard O2 DIN 9/16-18UNF
with 1 DIN quick-action coupling
++ ANTIMAGNETIC ++

Compact pressure regulator with contents gauge,
stepwise adjustable flow rate,
DIN quick-action coupling at a pressure of 450 kPa.

Gas type: Oxygen / O2

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 9

Outlet: 9/16-18 UNF (downwards, with tubing nipple)

1 quick-action coupling according to DIN 13260 Part 2 (upwards)

Flow rate: 0...15 l/min;

450 kPa / 120 l/min quick-action coupling

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type

O2-Oxygen
