

Hose, O2, 8m, DIN - NIST, Colour-Neutral



Manufacturer: Greggersen

Product 500165

Number:

Description

Low-Pressure Hose

according to DIN EN ISO 5359

Angled Probe with reverse flow limiter,
hose and NIST equipment connection

Specification:

Gas type: Oxygen / O2

Inlet: Probe according to DIN 13260 Part 2

Hose length: 8 m

Colour coding: Colour-Neutral

Outlet: gas-specific NIST screw connection

Properties

Gas Type O2-Oxygen