


Forano, Front Part, BS, MA, ISO



Manufacturer: Greggersen

Product 905161

Number:

Description

Terminal Unit FORANO - Front Part BS

Gas type: Compressed Air / MA

Inscription: medical air / MA

Gas-specific Socket Assembly for mounting
on Base Part FORANO.

Socket Assembly according to BS 5682.

Release via Release Ring.

Colour coding:

Socket Assembly coloured according to DIN EN ISO 9170-1 (ISO32).

Release Ring inscription: printed

Cartridge incl. sealing elements

Technical Data:

Gas type: Compressed Air / MA

Inscription: medical air / MA

Colour coding: according to DIN EN ISO 9170-1 (ISO32)

Front diameter: Ø 42mm

Scope of Delivery:

- gas-specific Socket Assembly
- Release Ring incl. printed inscription
- Cartridge with sealing elements

Properties

Gas Type	AIR Compressed Air
----------	--------------------
