

Flowmeter, DIN, AIR, 0-15l/min



Manufacturer: Greggersen

Product 900893

Number:

Description

Flowmeter, AIR, 0...15 l/min,

- PLUG-IN UNIT -

Medical flowmeter, pressure-compensated according to
DIN EN ISO 15002, with extension, fine
adjustment valve for metering the flow,
for coupling into a terminal unit

Gas type: Compressed Air - AIR

Connection: Probe according to DIN 13260 Part 2

Total length: 86.5 mm (probe front part to centre of flow outlet)

Outlet: 9/16-18 UNF

Flow rate: 0...15 l/min continuously adjustable

Inlet pressure: 450 kPa +/- 50 kPa

Properties

Gas Type	AIR Compressed Air
----------	--------------------
