

Flowmeter, NF, O₂, 0-15l/min



Manufacturer: Greggersen

Product 901500

Number:

Description

Flowmeter, O₂, 0...15 l/min, NF

- PLUG-IN UNIT -

Medical flowmeter, pressure-compensated according to DIN EN ISO 15002, with fine adjustment valve for metering the flow, for coupling into a terminal unit

Technical Data:

Gas type: Oxygen

Connection: Probe type NORME FRANCAISE

Outlet: 9/16-18 UNF

Flow rate: 0...15 l/min

Inlet pressure: 450 kPa +/- 50 kPa

Properties

Gas Type	O2-Oxygen
----------	-----------
