


Falke without Flow easy O2 DIN 9/16



Manufacturer: Greggersen

Product 910373

Number:

Description

Falke easy "RD", 450 kPa, 9/16-18UNF,

Compact pressure regulator with high-pressure cylinder connection to DIN 477, contents gauge, free flow at a pressure of 450 kPa.

Gas type: Oxygen / O₂

Connection: High-pressure cylinder connection to DIN 477, No. 9

Outlet: 9/16 -18UNF (Tubing Nipple) bottom

Output: 450 kPa

Properties

Gas Type O₂-Oxygen