

Forano, Base Part, UP/HW, DIN, BS, VAC



Manufacturer: Greggersen

Product 903252

Number:

Description

Terminal Unit Forano Base Part (Base Part) DIN/AGA/NF/BS

Gas-specific base part for installation in solid masonry walls (flush mounting) or in cavity walls.

Installation housing with threaded ring for subsequent height adjustment (up to 25 mm). Base part with brazed copper pipe; integrated non-return and maintenance valve.

Scope of delivery:

- gas-specific base part with Cu pipe
- installation housing with height adjustment
- chimney (for protection of the Cu pipe)
- dust protection caps

Technical data:

Gas type: Vacuum / VAC

Standard AGA: UTL (AGSS)

Installation type: Flush / cavity wall mounting

Inlet: Cu pipe 8x1 mm straight, top

Properties

Gas Type	VAC-Vacuum
----------	------------
