

Compressed-Air Vacuum Regulator Specht -90 DIN

Probe



Manufacturer: Greggersen

Product 900960

Number:

Description

Compressed-Air Vacuum Regulator Specht-90, Plug-in Unit DIN

Medical compressed-air vacuum regulator for vacuum generation, with vacuum gauge, fine adjustment and quick-shut-off valve, silencer and filter holder.

Technical Data:

Gas Type: Compressed Air / AIR

Inlet: Probe to DIN 13260 Part 2

Control Range: 0 to -85 kPa (0 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum hose barb

Properties

Gas Type AIR Compressed Air