

Kolibri 2 Probe DIN AIR Standard



Manufacturer: Greggersen

Product 906410

Number:

Description

KOLIBRI 2, AIR, 0...15 l/min

- Standard DIN PLUG-IN UNIT -

Flow Settings:

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selector mechanism

Medical compact flowmeter for connection
to a terminal unit, with handwheel - incrementally
adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised setting device and readability
of the set flow setting possible from various viewing directions
and angles.

Frontal operating concept -> enhanced safety when
setting the correct flow rate

Gas type: med. Compressed Air / AIR

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow Rate: 0...15 l/min at 450 kPa (4.5 bar) inlet pressure

Accessories required for operation:

900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	AIR Compressed Air
Inlet Standard	DIN
Variant	Standard