

Automatic Manifold MediControl MC2050E



Manufacturer: Greggersen

Product 326050

Number:

Description

Automatic Manifold MediControl 2050E

Flow Rate: 50 m³/h

Inlet pressure: max. 200 bar

Outlet pressure: 5 bar +/- 0.5 bar

Electr. connection: 230 V AC

System language: German

Mode: 2 cylinder banks

Gas type: ++++++

Consisting of:

1 x MC2050

1 x Pipe Set MC2050

Microprocessor-controlled gas supply system for 2-sided cylinder supply, design according to DIN EN ISO 7396-1, LCD display with operating status indication and plain text service messages, transmission of

messages via potential-free contacts, sensor-based pressure monitoring, 2-stage pressure reduction, redundant 2nd stage, pneum. priority switching in case of power failure, with pressure gauges and emergency supply option (NIST), fully assembled on mounting plate incl. cover