

Forano, Base Part, UP/HW, DIN, AIR



Manufacturer: Greggersen

Product 903251

Number:

Description

Terminal Unit Forano Base Part DIN/AGA/NF

Gas-specific base part for installation in solid masonry walls (flush-mounted) or in cavity walls.

Installation housing with threaded ring for subsequent height adjustment (up to 25 mm). Base part with brazed copper pipe; integrated non-return and maintenance valve.

Scope of delivery:

- gas-specific base part with Cu pipe
- installation housing with height adjustment
- protective sleeve (to protect the Cu pipe)
- dust protection caps

Technical data:

Gas type: Compressed air / AIR

Installation type: Flush-mounted / cavity wall mounting

Inlet: Cu pipe 8x1 mm straight, top

Properties

Gas Type	AIR Compressed Air
----------	--------------------
