

Vacuum Regulator Specht -60, AIR DIN Probe



Manufacturer: Greggersen

Product 901590

Number:

Description

Compressed-Air Vacuum Regulator Specht-60, DIN Plug-in Unit

Medical compressed-air vacuum regulator for vacuum generation, with vacuum gauge, fine adjustment and quick-action shut-off valve, silencer and filter mount.

Technical Data:

Gas type: Compressed air / AIR

Inlet: Probe according to DIN 13260 Part 2

Control range: 0 to -60 kPa (0 to -0.60 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air
-----------------	--------------------