

Portable Emergency Supply Habicht Kolibri Comfort



Manufacturer: Greggersen

Product 902698

Number:

Description

Emergency Supply Unit,
portable, without oxygen cylinder
- with KOLIBRI Compact Flowmeter -

Emergency supply unit, oxygen flowmeter
and ejector for vacuum generation
consisting of:

Stainless steel carrying frame with stand, to accommodate
a 2 L or 3 L O₂ cylinder, suitable
for attaching to medical rail or bed,
central supply block with NIST inlet,
consisting of:
KOLIBRI compact flowmeter 0...15 l/min with flow settings
for metering the flow rate and ejector with vacuum gauge
and fine adjustment valve;

Gas type: Oxygen / O₂

Connection: NIST according to DIN EN ISO 18082

Outlet: Tubing Nipple / vacuum nipple

Performance: FLOW : 0...15 l/min

VAC: adjustable between

0 and -85 kPa (0 to -0.85 bar)