

Flowmeter, DIN, O₂, 0-15l/min extended



Manufacturer: Greggersen

Product 900754

Number:

Description

Flowmeter, O₂, 0...15 l/min

with extension

PLUG-IN UNIT

Medical flowmeter, pressure-compensated according to
DIN EN ISO 15002, with extension, fine
adjustment valve for metering the flow,
for connecting to a terminal unit

Gas type: Oxygen - O₂

Connection: Probe according to DIN 13260 Part 2

Total length: 131.50 mm (probe front part to centre of flow outlet)

Outlet: 9/16 -18 UNF

Flow rate: 0...15 l/min continuously adjustable

Inlet pressure: 450 kPa +/- 50 kPa

Properties

Gas Type	O2-Oxygen
----------	-----------
