

Kolibri 2 Double Probe DIN O2 Standard



Manufacturer: Greggersen

Product 906408

Number:

Description

KOLIBRI 2, O₂, 0...15 l/min

- Standard - Double - DIN PLUG-IN UNIT -

Flow settings (2x):

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Double Compact Flowmeter KOLIBRI 2

with fixed orifices and selector device

Medical compact flowmeter for coupling
into a terminal unit, with handwheel - incrementally
adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and readability
of the set flow setting from various viewing directions
and angles.

Frontal operating concept -> enhanced safety when
setting the correct flow rate

Gas type: Oxygen / O₂

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow rate: 0...15 l/min at 450 kPa (4.5 bar) inlet pressure

Accessories required for operation:

2x 900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	O ₂ -Oxygen
Inlet Standard	DIN
Variant	Standard Double