


Falke with Flow Neo O2 DIN 9/16



Manufacturer: Greggersen

Product 910353

Number:

Description

Falke with Flow, NEO

Compact pressure regulator contains gauge;
stepwise adjustable flow rate.

Gas type: Oxygen / O₂

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 9

Outlet: 9/16-18 UNF (tubing nipple) bottom

Flow rate: 0...1 l/min

Flow settings: 0/ 0.1/ 0.15/ 0.2/ 0.25/ 0.3/ 0.4/ 0.5/ 0.6/ 0.8/ 1.0 l/min

Properties

Gas Type

O₂-Oxygen