


Falke with Flow Paediatric O2 DIN Tubing Nipple



Manufacturer: Greggersen

Product 910202

Number:

Description

Falke with Flow, Paediatric,
Compact pressure regulator with , contents gauge;
incrementally adjustable flow rate.

Gas type: Oxygen / O₂

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 9

Outlet: Tubing Nipple bottom

Flow Rate: 0...5 l/min.

Flow Settings: 0/ 0.1/ 0.3/ 0.5/ 0.6/ 0.8/ 1/ 2/ 3/ 4/ 5 l/min

Properties

Gas Type O₂-Oxygen