

Kolibri 2 DIN Probe O2 Standard+



Manufacturer: Greggersen

Product 906401

Number:

Description

KOLIBRI 2, O₂, 0...30 l/min

- Standard+ DIN PLUG-IN UNIT -

Flow Settings:

0 / 2 / 3 / 4 / 6 / 8 / 12 / 15 / 20 / 25 / 30 l/min

Compact Flowmeter KOLIBRI 2

with fixed orifices and selection device

Medical compact flowmeter for coupling
into a terminal unit, with handwheel – incrementally
adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and readability
of the set flow setting possible from various viewing directions
and angles.

Frontal operating concept -> greater safety when
setting the correct flow rate

Gas Type: Oxygen / O₂

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow Rate: 0...30 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:

900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	O ₂ -Oxygen
Inlet Standard	DIN
Variant	Standard+