


## Pressure Regulator Unit, AIR, Size 2, 5bar



**Manufacturer:** Greggersen

**Product** 325815

**Number:**

### Description

Pressure regulator unit for reducing the receiver pressure to 5 bar, design based on DIN EN ISO 7396-1, each pressure regulator stage duplicated with shut-off valves, safety valve, contact gauge, fully mounted on a mounting plate.

Size: 2

Capacity: 300 Nm<sup>3</sup>/h

Gas type: Compressed air / AIR

Pressure: 500 kPa (5 bar)

Contact gauge: 400 / 600 kPa (4 bar / 6 bar)

### Properties

**Gas Type** AIR Compressed Air