


Kolibri, Rail O2 Standard+ - NIST



Manufacturer: Greggersen

Product 902601

Number:

Description

Kolibri, O₂, 0...30 l/min

RAIL-MOUNTED UNIT NIST

Medical compact flowmeter

according to DIN EN ISO 15002; flow rate is set via
handwheel in 10 flow settings;

gas inlet NIST integrated into rail clamp;

outlet 9/16-18 UNF and tubing nipple

Gas type: Oxygen O₂

Connection: NIST according to DIN EN ISO 18082

Outlet: 9/16-18 UNF / Tubing Nipple

Flow rate: 0...30 l/min

Flow settings: 0 / 0.5 / 1 / 2 / 4 / 6 / 8 / 10 / 12 / 15 / 30+ l/min

Properties

Gas Type	O2-Oxygen
----------	-----------
