

## Falke w. Flow Stand.+ O2 DIN 9/16 +1DIN



**Manufacturer:** Greggersen

**Product** 910356

**Number:**

### Description

Falke with flow, Standard+ with 1 quick-release coupling,  
Compact pressure regulator with contents gauge;  
step-adjustable flow rate,  
DIN quick-release coupling at a pressure of 450 kPa.

Gas type: Oxygen / O2

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 9

Outlet: 9/16-18 UNF (Tubing Nipple) / 1 quick-release coupling according to DIN 13260-2

Flow rate: 0...30 l/min 9/16-18 UNF; 120 l/min quick-release coupling

Flow settings: 0 / 0.5/ 1/ 2/ 4/ 6/ 8/ 10/ 12/ 15/ 30+ l/min

### Properties

**Gas Type**

O2-Oxygen