

Vacuum Regulator Spatz -16 DIN Probe



Manufacturer: Greggersen

Product 900899

Number:

Description

Vacuum Regulator Spatz-16, Plug-in Unit DIN

Medical vacuum regulator with vacuum gauge,
fine adjustment and quick-acting shut-off valve.

Technical data:

Gas type: Vacuum / VAC

Inlet: Probe according to DIN 13260 Part 2

Control range: 0 to -16 kPa (0 to -160 mbar)

Outlet: 9/16-18 UNF with vacuum tubing nipple