


Falke without Flow N2O DIN +1DIN (small cylinder)



Manufacturer: Greggersen

Product 910384

Number:

Description

Falke, without flow with 1 quick-release coupling,
Compact pressure regulator contains gauge;
with one quick-release coupling at a pressure of 450 kPa

Gas type: Nitrous oxide / N₂O

Connection: Cylinder connection for gas cylinders according to DIN 477, No. 12 (small cylinder)

Outlet: Quick-release coupling according to DIN 13260 Part 2, bottom

Output: 450 kPa

Properties

Gas Type N₂O-Nitrous Oxide