


Ventus 22 mm Valve Insert N2O



Manufacturer: Greggersen

Product 900783

Number:

Description

Valve insert for Ventus 22 mm

Gas type : N2O

Design complies with DIN EN ISO 7396-1

Valve assembly consisting of:

- Shut-off valve DN 20, G3/4"
- Cu pipe 22mm inlet / outlet
- Pressure gauge 0...16 bar
- Pressure sensor -1...6
- NIST emergency inlet point
- Physical disconnection point
- Gas type labelling

Properties

Fulfilled Standard	DIN EN ISO 7396-1
Gas Type	N2O-Nitrous Oxide
NIST Emergency Inlet Point	1 NIST Emergency Inlet Point