


Falke without Flow N2O DIN +1DIN (large cylinder)



Manufacturer: Greggersen

Product 910378

Number:

Description

Falke without Flow with a DIN Quick Coupling

Compact pressure regulator with contents pressure gauge;

DIN quick coupling at a pressure of 450 kPa.

Gas type: Nitrous Oxide / N2O

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 11 (large cylinder)

Outlet: Quick Coupling DIN 13260 Part 2, top

Output: 450 kPa

Properties

Gas Type N2O-Nitrous Oxide