


Hose, AIR, 5m, DIN - DIN Coupling Colour-Neutral



Manufacturer: Greggersen

Product 902531

Number:

Description

Low-Pressure Hose
to DIN EN ISO 5359

Angled Probe with reverse flow limiter,
hose and coupling to DIN

Specification:

Gas type: Compressed Air / AIR

Inlet: Probe to DIN 13260 Part 2

Hose length: 5.0 m

Colour coding: Colour-Neutral

Outlet: Coupling to DIN 13260 Part 2

Properties

Gas Type AIR Compressed Air