

Kolibri 2 Double Probe DIN AIR Standard



Manufacturer: Greggersen

Product 906414

Number:

Description

KOLIBRI 2, AIR, 0...15 l/min

- Standard - Double - DIN PLUG-IN UNIT -

Flow settings (2x):

0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Double Compact Flowmeter KOLIBRI 2

with fixed orifices and selector mechanism

Medical compact flowmeter for connection

to a terminal unit, with handwheel - incrementally

adjustable flow rate, outlet 9/16-18 UNF thread

Ergonomically optimised adjustment device and readability

of the selected flow setting from various viewing directions

and angles.

Frontal operating concept -> greater safety when
setting the correct flow rate

Gas type: Med. Compressed Air / AIR

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow rate: 0...15 l/min at inlet pressure 450 kPa (4.5 bar)

Accessories required for operation:

2 x 900619 - Hose connection 9/16-18 UNF

Properties

Gas Type	AIR Compressed Air
Inlet Standard	DIN
Variant	Standard Double