

Flowmeter, DIN, CO2, 0-12l/min



Manufacturer: Greggersen

Product 500011

Number:

Description

Flowmeter, Standard, CO2

Plug-in Unit

Medical flowmeter, pressure-compensated according to
DIN EN ISO 15002, with fine adjustment valve for
flow rate dosing, for coupling into
a terminal unit

Technical Data:

Gas type: Carbon dioxide / CO2

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow rate: 0...12 l/min

Properties

Gas Type CO2-Carbon Dioxide